

## HALO H<sub>2</sub>O Trace Level Moisture Analyzer for Nitrous Oxide

 GASES & CHEMICALS
 CEMS
 ENERGY
 ATMOSPHERIC
 SEMI & HB LED
 SYNGAS
 LABORATORY

## Designed for trace level moisture analysis in nitrous oxide, the HALO H<sub>2</sub>O offers:

- Low parts per billion (ppb) moisture detection capability in medical N<sub>2</sub>O
- Absolute measurement (freedom from calibration gases)
- Wide dynamic range over four orders of magnitude
- Low cost of ownership and operational simplicity
- Clean technology no external calibration gases required
- Compact analyzer footprint

The HALO  $H_2O$  analyzer provides users with the unmatched accuracy, reliability, speed of response and ease of operation that users of Tiger Optics analyzers have come to expect. Featuring Tiger Optics' powerful Cavity Ring-Down Spectroscopy-based moisture sensor in a compact and economical package, this versatile analyzer allows users to measure trace moisture in nitrous oxide (N<sub>2</sub>O) and other gases. Users can also expect significant cost savings, with no need for periodic sensor maintenance, span calibrations, purifier replacements, and pump rebuilds. As a result, the HALO  $H_2O$  analyzer is ideally suited to many applications where moisture measurement is extremely critical. These applications include process control, continuous quality control, and cylinder analysis for medical and calibration gases and more.



## **HALO H<sub>2</sub>O** Trace Level Moisture Analyzer for Nitrous Oxide



Performance		
Operating range	See table below	
Detection limit (LDL,	See table below	
24 h peak-to-peak variation)		
Sensitivity (3o)	See table below	
Precision (1 $\sigma$ , greater of)	± 0.75% or 1/3 of Sensitivity	
Accuracy (greater of)	± 4% or 1/2 of LDL	
Speed of response	< 1 minute to 90%	
Environmental conditions	10°C – 40°C	
	30% – 80% RH (non-condensing)	
Storage temperature	-10°C – 50°C	

## **Gas Handling System and Conditions**

Wetted materials	316L stainless steel		
	10 Ra surface finish		
Gas connections	1/4" male VCR inlet and outlet		
Leak tested to	1 x 10 <sup>-9</sup> mbar l / sec		
Inlet pressure	10 – 125 psig (1.7 – 9.6 bara)		
Flow rate	Up to 1.8 slpm		
Sample gases	Most inert, toxic, passive		
	and corrosive matrices		
Gas temperature	Up to 60°C		

Dimensions	H x W x D [in (mm)]		
Standard sensor	8.75 x 8.5 x 23.6 (222 x 216 x 599)		
Sensor rack	8.75 x 19 x 23.6 (222 x 483 x 599)		
(fits up to two sensors)			
Weight			
Standard sensor	28 lbs (12.7 kg)		
Electrical			
Alarm indicators	2 user programmable		
	1 system fault		
	Form C relays		
Power requirements	90 – 240 VAC, 50/60 Hz		
Power consumption	40 Watts max.		
Signal output	Isolated 4-20 mA per sensor		
User interfaces	5.7" LCD touchscreen		
	10/100 Base-T Ethernet		
	802.11g Wireless (optional)		
	RS-232		

Performance, H <sub>2</sub> O:	Range	LDL	Sensitivity
In Nitrogen	0 – 20 ppm	3 ppb	2.2 ppb
In Nitrous Oxide (N <sub>2</sub> O)	0 – 20 ppm	10 ppb	7.5 ppb

Contact us for additional analytes and matrices. U.S. Patent # 7,277,177



